

# CONVENTIONAL HIGH PERFORMANCE LIPOSUCTION UNIT *MODEL VAC 1288-R3*



- \* EXCELLENT VACUUM PRESSURE
  - \* ULTRA-SILENT !!!.
- \* VACUUM POWER 1/3 HP.
- \* CONTINUOUS SUCTION.
- \* BIOLOGIC CONTROL CANISTER.
- \* IDEAL LIPO-SCULPTURE SYSTEM WITH OUR VIBRO-LIPO-SET-R3



## TECHNICAL FEATURES

- ✓ Conventional high performance liposuction unit.
- ✓ Ultra-silent.
- ✓ Suction capacity 3.000 cc. with biologic control canister.
- ✓ Maximum vacuum pressure: 1 Atm., at see level.
- ✓ Suction power 1/3 HP.
- ✓ Vacuum pressure indicated at the panel gauge, with adjusting control.
- ✓ Easy to transport with four wheels system.
- ✓ Possibility of using as tower configurating our Vibro-Lipo-Set-R3
- ✓ Voltage 110 / 120 Volts A.C. other voltages upon request.