



INFILTRATION PUMP MODEL IFP 2000-D



- * SUBCUTANEOUS INFILTRATION SYSTEM THROUGH A DIGITAL ELECTRONIC CONTROLLED THREE ROLLER PUMP.
- * MOST UPDATED LIPO-SCULPTURE HIGH FLOW INFILTRATION SYSTEM.
- *TWO OPERATION MODES : STANDARD AND VOLUMETRIC PROGRAMMABLE.
- * EASES PROCEDURE MINIMIZING SURGICAL TIME AND REDUCING COSTS.
- * IDEAL LIPO-SCULPTURE SYSTEM WITH OUR VIBRO-LIPO-SET-PC AND VIBRO-LIPO-SET-R3

TECHNICAL FEATURES

- ✓Two operation modes: Standard and Volumetric Programmable.
- ✓Infiltration rate (flow) and volume to infiltrate adjustable through front panel membrane keyboard and Visualized at alphanumerical digital display.
- ✓Infiltration rate (flow) adjustment range: 400 cc / min. – 800 cc / min.
- ✓Volume to infiltrate (Volumetric mode) programming range: 50 cc to 1000 cc in 50 cc steps.
- ✓Total infiltrated volume computing and showed at front panel alphanumerical display.
- ✓Program for easy insertion and remove of pumping silastic hose segment.
- ✓Stored activation through pedal and / or manual through the front panel membrane keyboard.
- ✓Lightweight and easy to transport.
- ✓Possibility of using on table or on trolley with our Vibro-Lipo-Set-PC and Vibro-Lipo-Set-R3
- ✓Voltage 110 / 120 Volts A.C. other voltages upon request.