



# PORTABLE HIGH PERFORMANCE LIPOSUCTION UNIT MODEL VAC 1288-L

- \* MULTIPURPOSE USE
- \* EXCELLENT VACUUM PRESSURE
- \* HIGH FLOW
- \* LIGHTWEIGHT AND EASY TO TRANSPORT
- \* IDEAL FOR ALL LARGE AND SMALL SURGICAL PROCEDURES
- \* CONTINUOUS SUCTION
- \* BIOLOGIC CONTROL CANISTER



## TECHNICAL FEATURES

- ✓ Maximum vacuum pressure: 1 Atm. at sea level.
- ✓ Frontal control for adjusting vacuum pressure.
- ✓ Suction capacity: Biologic control canister of 1500 CC., 3000 CC. upon request.
- ✓ High flow Continuous suction.
- ✓ Voltage 110 / 120 Volts. A.C. other voltages upon request.