



DIGITAL SHAVER ARTHRO-DRIVE CHONDROTOME MODEL AH 1507-D

- * IDEAL EQUIPMENT FOR ARTHROSCOPY
- * EXCELLENT SCRAPING QUALITY
- * ERGONOMIC AND LIGHTWEIGHT HANDPIECE
- * LIQUID CRYSTAL ALPHANUMERICAL DISPLAY



TECHNICAL FEATURES:

- ✓ Adjustable speed range between 500 and 4000 RPM visualized at the front panel.
- ✓ Liquid crystal alphanumerical display which shows motion direction: left, right, oscillatory , and the actual spin speed in RPM.
- ✓ Ergonomic lightweight hand-piece with suction electronic control through a membrane knob.
- ✓ Hand piece engine speed and spin direction electronic control through a double pedal, with alarm function to warn about overloading or damage of hand-piece cable.
- ✓ Stryker blades compatible.
- ✓ Voltage 110 / 120 Volts A.C. other voltages upon request.